

# Android MacroTelectLink Development Guide

---

## Introduction

This guide will teach you how to use MacroTelectLink SDK to write Android applications that can acquire brainwave data from MacroTelect 's Hardware (BrainLink Pro & BrainLink Lite) . This will enable your Android apps to receive and use brainwave data such as BrainWave and Gravity acquired via Bluetooth, MacroTelect 's Hardware and File source are encapsulated as MacroTelectLink SDK for Android supports upgrading Hardware

**Function:** Receive brainwave data.

### Supported Device:

- Bluetooth 4.0 BLE
  - BrainLink\_Pro
- Bluetooth 3.0
  - BrainLink\_Lite
  - Mind Link

**Supported Android Version:** - Android 4.3 +

## Your First Project: MacroTelectLinkDemo (Android studio)

1. Copy the SDK to the project's libs folder and add dependencies in build.gradle.

```
dependencies{  
    ..  
    implementation files('libs/MacroTelectLink_V1.4.3.jar')  
}
```

2. Add permissions in AndroidManifest.xml ``

```
**3.Receive brainwave data from MacroTelect SDK.**  
```java  
  
//check location permissions  
if((checkSelfPermission(Manifest.permission.ACCESS_COARSE_LOCATION)!=PackageManager.PERMISSION_GRANTED))  
    //request permissions  
    requestPermissions(new String[]{Manifest.permission.ACCESS_COARSE_LOCATION,Manifest.permission.ACCESS_FINE_LOCATION},1)
```

```

bluemanage = LinkManager.init(this);
bluemanage.setDebug(true); //whether to print

//set brainwave data callback
bluemanage.setMultiEEGPowerDataListener(new EEGPowerDataListener() {
    @Override
    public void onBrainWavedata(String mac, BrainWave brainWave) {
        Log.e(mac, brainWave.toString()); //receive brainwave data
    }
    @Override
    public void onRawData(String mac, int raw) {
        //receive raw data
    }
    @Override
    public void onGravity(String mac, Gravity gravity) {
        //receive gravity data
    }
    @Override
    public void onRR(String mac, ArrayList<Integer> rr, int oxygen) {
        //Receiving data: RR intervals and blood oxygen percentage
    }
});

//connection status recall
bluemanage.setOnConnectListener(new OnConnectListener() {
    @Override
    public void onConnectionLost(BlueConnectDevice blueConnectDevice) {
        // lost connected Bluetooth
    }
    @Override
    public void onConnectStart(BlueConnectDevice blueConnectDevice) {
        //trying to connect Bluetooth
    }
    @Override
    public void onConnectting(BlueConnectDevice blueConnectDevice) {
        //connecting...
    }
    @Override
    public void onConnectFailed(BlueConnectDevice blueConnectDevice) {
        //Bluetooth disconnected
    }

    @Override
    public void onConnectSuccess(BlueConnectDevice blueConnectDevice) {
        String mac = blueConnectDevice.getAddress();
        String connectType = blueConnectDevice.isBleConnect ? " 4.0 " : " 3.0 ";
        Log.e(TAG, "connected device name:" + blueConnectDevice.getName() + " mac: " +
            mac);
    }
});

```

```

    }

    @Override
    public void onError(Exception e) {
        Log.e(TAG, "connect Bluetooth error");
        e.printStackTrace();
    }
});

bluemanage.setMaxConnectSize(1);//Set the maximum of connections
bluemanage.setConnectType(LinkManager.ConnectType.ALLDEVICE);//set connect type
bluemanage.setWhiteList("BrainLink_pro,BrainLink_Lite");//set whitelist
bluemanage.startScan();//start scan and connect

```

## MacrotellectLinkAPI Reference

---

### Brainwave Reference

---

Brainwave information entity class, a real-time analysis of the chip to the software that obtain detailed brainwave information and power data.

**Note:** \ Connect to BrainLink\_Lite and Mind Link, which only have : signal, att, med, delta, theta, lowAlpha, highAlpha, lowBeta, highBeta, lowGamma, middleGamma, heartRate and temperature.

#### Basic Brainwave Data

- **Signal** (It represents the signal value of the Macrotellect's hardware. When the signal is 0, it means that the hardware has been put on, and when the signal is 200, it means that hardware is connected to the phone via Bluetooth.)
- **att**(Attention)
- **med**(Relaxation)
- **delta**
- **theta**
- **lowAlpha**,
- **highAlpha**
- **lowBeta**
- **highBeta**
- **lowGamma**
- **middleGamma**
- **ap**(Appreciation)
- **batteryCapacity**(BatteryCapacity)
- **heartRate**
- **temperature**

# GravityReference

---

Gravity information entity class, a real-time analysis of the chip to software that obtain detailed gravity information data.

**Note:** Connect to BrainLink\_Lite has no function of receiving gravity data.

**Gravity Data** - X value: gravity value in The x axis(Pitching angle) - Yvalue:gravity value in The yaxi(Yaw angle) - Z value: gravity value in The z axis(Roll angle)

## EEGPowerDataListener Reference

---

The class for receiving real-time parse chip data needs to be implemented by the user, which receive brainwave、gravity and raw EEG data. Set the data parsing result callback in the `setEegPowerDataListener (OnConnectListenerOnConnectListener)` of the Linkmanager class.

### Method

**void onBrainWavedata(String mac,BrainWavebrainWave);**

Receive parsed brainwave data. - mac: mac address of Brainwave device - brainWave:brainwave data

**void onGravity(String mac,Gravity gravity);**

Receive parsed gravity data.

- mac: mac address of Brainwave device
- gravity:gravity data

**void onRawData(String mac,int raw);**

Receive raw EEG data.

- mac: mac address of Brainwave device
- raw:raw EEG data

**void onRR(String mac, ArrayList rr, int oxygen);**

Receiving data: RR interval and blood oxygen percentage.

- mac: mac address for BrainLink device
- rr: RR intervals
- oxygen: blood oxygen percentage

## OnConnectListenerReference

---

This class is the connection state monitoring interface between Bluetooth and brainwavedevice needs to be implemented by the user. Set up the connection listener in the Linkmanager's setOnConnectListener (OnConnectListener).

## Method

**void onConnectStart(BlueConnectDevice blueConnectDevice);**

trying to connect

**void onConnectting(BlueConnectDevice blueConnectDevice);**

connecting

**void onConnectFailed(BlueConnectDevice blueConnectDevice);**

disconnected

**void onConnectSuccess(BlueConnectDevice blueConnectDevice);**

connected

**void onConnectionLost(BlueConnectDevice blueConnectDevice);**

Lost connection (disconnected from connected state) **void onError(Exception e);**

connection error

# LinkManagerReference

---

This class handles the interaction between the MacroTelect's hardware and Bluetooth devices.

## Method

**public static LinkManager init(Context context)**

Initialization (singleton)

**public void setDebug(boolean isDebug)**

whether to print log (no print setting by default)

**public void setMaxConnectSize(int count)**

Set the maximum number of connections (1 is set by default)

**public void setConnectType(ConnectType connectType);**

set the type of connection

- **ConnectType. ONLYCLASSBLUE** Only connect devices by class Bluetooth, you need to manually pair them first.
- **ConnectType. ONLYBLEBLUE** Only connect devices by BLE Bluetooth
- **ConnectType. ALLDEVICE** Allow both ways to connect devices

**getConnectSize();**

receive the number of connected devices

**public void setWhiteList(String whiteList)**

Set whitelist, only allow to connect whitelist. Please use ' , ' to separate names of connected multiple devices.

**public void setOnConnectListener(OnConnectListener onConnectListener)**

Set the Bluetooth connection status callback

**public void setEegPowerDataListener(EEGPowerDataListener eegPowerDataListener)**

Set the brainwave data receiving callback

## **Change Records**

**Add blood oxygen percentage**